

Arterial Vascular Impedance Calculation

AIM: To calculate Arterial Vascular Impedance, $Z(j\Omega) = \frac{P(j\Omega)}{Q(j\Omega)}$, where P is intra-arterial pressure (in mmHg) and Q is blood volume-flow (in ml/s).

Currently, the data (1,2,3) can be collected simultaneously using CMCdaq without any time- synchronization issue.

STEP1: The following are the data collected simultaneously:

1. Intra-arterial pressure ($p[n]$) at 4 KHz.

2. Video recorded at 30fps.

3. Audio recorded at 4KHz



From Ultrasound Machine

STEP2: Conversion to frequency domain:

2.1. $p[n] \rightarrow \boxed{\text{DFT}} \rightarrow P(j\Omega)$, N= 4000 samples.

2.2. Calculation of volume-flow: $q[n]$

2.2.1. Frame-wise area calculation: $a[n]$, $n = \text{frame no.}$

2.2.2. For one frame duration, no: of audio samples = $4000/30 = 133$ samples: $x[n]$, $N=133$

2.2.3. $x[n] \rightarrow \boxed{\text{DFT}} \rightarrow X(j\Omega) \rightarrow \boxed{\text{Averaging}} = \text{average doppler shifted freq } (f_D) \rightarrow \text{conversion to average velocity, } v[n]$.

2.2.4. Hence, calculation of volume flow, $q[n] = a[n] \cdot v[n]$

2.3. Conversion of $q[n]$ to frequency domain:

$q[n] \rightarrow \boxed{\text{DFT}} \rightarrow Q(j\Omega) \rightarrow \boxed{\text{Zero padding to 4000 samples}} \rightarrow \boxed{\text{IDFT}} \rightarrow q[n]$, $N=4000$

2.4. Arterial Vascular Impedance, $Z(j\Omega) = \frac{P(j\Omega)}{Q(j\Omega)}$