

Arterial Vascular Impedance Calculation

AIM: To calculate Arterial Vascular Impedance, $Z(j\Omega) = \frac{P(j\Omega)}{Q(j\Omega)}$, where P is intra-arterial pressure (in mmHg) and Q is blood volume-flow (in ml/s).

Currently, the data (1,2,3) can be collected simultaneously using CMCdaq without any time- synchronization issue.

STEP1: The following are the data collected simultaneously:

1. Intra-arterial pressure ($p[n]$) at 4 KHz.

2. Video recorded at 30fps.

3. Audio recorded at 4KHz



From Ultrasound Machine

STEP2: Conversion to frequency domain:

2.1. $p[n] \rightarrow$ **DFT** $\rightarrow P(j\Omega)$, $N= 4000$ samples/s.

2.2. Calculation of volume-flow: $q[n]$

2.2.1. Frame-wise area calculation: $a[n]$, $n= \text{frame no.}$

2.2.2. For one frame duration, no: of audio samples = $4000/30 = 133$ samples/frame: $x[n]$, $N=133$

2.2.3. $x[n] \rightarrow$ **DFT** $\rightarrow X(j\Omega) \rightarrow$ **Averaging** = average doppler shifted freq (f_D) \rightarrow conversion to average velocity, $v[n]$,

using, $f_D = \frac{2 \cdot f \cdot v \cos \theta}{c}$, f is the probe txn freq=11MHz, c is velocity of sound in blood=1580 m/s, $\theta = 60^\circ$.




2.2.4. Hence, calculation of volume flow, $q[n] = a[n] \cdot v[n]$

2.3. Conversion of $q[n]$ to frequency domain:

$q[n] \rightarrow$ **DFT** $\rightarrow Q(j\Omega) \rightarrow$ **Zero padding to 4000 samples** \rightarrow **IDFT** $\rightarrow q[n]$, $N=4000$

2.4. Arterial Vascular Impedance, $Z(j\Omega) = \frac{P(j\Omega)}{Q(j\Omega)}$

Summary of Calculations from Video (10fps) and audio (4KHz) collected simultaneously:



NB: Video was recorded at 10fps in CMCdaq. We need to standardize data-collection procedure (for stable angle & video).

